

ACCOUNT INFO

Name: _____ Account #: _____ PO# / MR#: _____

Shipping address: _____ City: _____ State: _____ Zip: _____

PATIENT INFO

Last name: _____ First name: _____

Age/DOB: _____ Weight (Required): _____ Height: _____ Male Female

Shoe type: Dress Casual Athletic Other: _____

Shoe size (Required): _____ Shoe width: _____ Diagnosis: _____

ORDER INFO

Make _____ pairs

Make from prior Rx: (# _____)

Ship to patient (please provide address)

Left only Right only

NORTHWEST SELECT ORTHOTICS™ Select one. Specs on reverse. **Patient weight required.**

Activity/use

- Athletic
 - Low profile shoes
 - Standard / high profile shoes
 - Standard / high profile shoes + lateral stability
- Everyday
- Dress (Heel below 1.5")
- Accommodative (Control + comfort)

Foot pain/pathology

- AAF/PTTD
- Heel spurs
- Sesamoiditis/hallux limitus
- Plantar fasciitis
- Amputee

YOU BUILD ORTHOTICS™ Patient weight required. Underlined options are standard/default. Rx forms with missing/incomplete info will automatically revert to standard/default options.

Shell/foundation

Polypropylene

High-quality polypropylene routed to NWPL's exacting specifications for excellent performance.

Semi-flexible Semi-rigid Rigid Extra rigid

Semi-flexible Semi-rigid Rigid

How it works

Select flex based on desired in-shoe feel & treatment goals

Ensure patient weight is provided

Shell crafted based on selected flex and patient information

Arch height (shell)

- Very high (Very intimate to arch)
- High (Intimate to arch)
- Standard
- Low (Less intimate to arch)

Top cover length and material

- 3/4
- To sulcus
- Heel to toe

Vinyl Charcoal Dark brown Gray-tan Blue-gray Leather (black)

Neoprene 1/16" (1.5mm) 1/8" (3.0mm) Microsuede Black Gray

Firm EVA 1/16" (1.5mm) 1/8" (3.0mm) Soft EVA 1/16" (1.5mm) 1/8" (3.0mm)

1/16" top cover requires extension

X-Guard (bottom cover)

- Standard
- Sulcus
- Full
- Complete

Cushion

- 3/4
- To sulcus
- Heel to toe

P-Foam 1/16" (1.5mm) 1/8" (3.0mm)

Soft EVA 1/16" (1.5mm) 1/8" (3.0mm)

Firm EVA 1/16" (1.5mm) 1/8" (3.0mm)

Extension

- Full
- Sulcus

P-Foam 1/16" (1.5mm) 1/8" (3.0mm)

Soft EVA 1/16" (1.5mm) 1/8" (3.0mm)

Firm EVA 1/16" (1.5mm) 1/8" (3.0mm)

DETAILS/SPECS FOR NORTHWEST SELECT ORTHOTICS™




































































Activity/use

Everyday	Semi-rigid shell / charcoal vinyl top cover (heel to toe) / 1/16" firm EVA extension (full) / X-Guard (standard)
Dress - Heel below 1.5"	Semi-rigid shell / black microsuede top cover (to sulcus) / 1/16" P-Foam extension (sulcus-length) / X-Guard (standard)
Athletic - Low profile shoes	Rigid shell / charcoal vinyl top cover (3/4-length) / X-Guard (standard)
Athletic - Standard / high profile shoes	Rigid shell / 1/16" neoprene top cover (heel to toe) / 1/16" firm EVA extension (full) / deep heel cup / X-Guard (standard)
Athletic - Standard / high profile shoes + stability	Rigid shell / 1/16" neoprene top cover (heel to toe) / 1/16" firm EVA extension (full) / deep heel cup / heel post (vertical) / X-Guard (standard)
Accommodative - Control + comfort	Semi-flexible shell / charcoal vinyl top cover (heel to toe) / 1/16" P-Foam cushion (heel to toe)

Foot pain/pathology

AAF/PTTD	Rigid shell (AAF/PTTD design) / charcoal vinyl top cover (3/4-length) / low arch height / medial heel skive / deep heel cup / integrated extrinsic forefoot post (3 degrees varus) / integrated heel post (vertical) / high medial overlay
Plantar fasciitis	Semi-rigid shell / charcoal vinyl top cover (heel to toe) / 1/16" firm EVA extension (full) / medial heel skive / deep heel cup / high arch height / X-Guard (standard)
Sesamoiditis/hallux limitus	Semi-rigid shell / charcoal vinyl top cover (heel to toe) / 1/16" Firm EVA extension (full) / Reverse Morton's extension
Heel spurs	Semi-rigid shell / charcoal vinyl top cover (heel to toe) / 1/16" P-Foam cushion (heel to toe) / 1/16" firm EVA extension (full) / deep heel cup / heel cushion
Amputee	Semi-rigid shell (amputee design) / 1/16" soft EVA top cover (heel to toe) / 1/16" firm EVA extension (full)

ORTHOTIC OPTIONS = Bilateral

<h3>Posting</h3> <p>Integrated heel post</p> <table style="width:100%; border: none;"> <tr> <td style="width:50%; border: none;"> <p>Left</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p> </td> <td style="width:50%; border: none;"> <p>Right</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p> </td> </tr> </table> <p style="text-align: center; font-size: small;">Maximum of 4°</p> <p style="font-size: x-small;">Integrated heel posts are crafted as part of the shell. Standard heel posts are available for an added cost. Request standard posts in notes.</p> <p>Heel lift (requires heel post)</p> <table style="width:100%; border: none;"> <tr> <td style="width:50%; border: none;"> <p>Left</p> <p>1/8" (3.0mm)</p> <p>2/8" (6.0mm)</p> <p>3/8" (9.0mm)</p> <p>4/8" (12.0mm)</p> </td> <td style="width:50%; border: none;"> <p>Right</p> <p>1/8" (3.0mm)</p> <p>2/8" (6.0mm)</p> <p>3/8" (9.0mm)</p> <p>4/8" (12.0mm)</p> </td> </tr> </table> <p>Integrated forefoot post</p> <table style="width:100%; border: none;"> <tr> <td style="width:50%; border: none;"> <p>Left</p> <p><input type="checkbox"/> Extrinsic</p> <p><input type="checkbox"/> Intrinsic</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p> </td> <td style="width:50%; border: none;"> <p>Right</p> <p><input type="checkbox"/> Extrinsic</p> <p><input type="checkbox"/> Intrinsic</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p> </td> </tr> </table> <p style="text-align: center; font-size: small;">Maximum of 6°</p>	<p>Left</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p>	<p>Right</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p>	<p>Left</p> <p>1/8" (3.0mm)</p> <p>2/8" (6.0mm)</p> <p>3/8" (9.0mm)</p> <p>4/8" (12.0mm)</p>	<p>Right</p> <p>1/8" (3.0mm)</p> <p>2/8" (6.0mm)</p> <p>3/8" (9.0mm)</p> <p>4/8" (12.0mm)</p>	<p>Left</p> <p><input type="checkbox"/> Extrinsic</p> <p><input type="checkbox"/> Intrinsic</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p>	<p>Right</p> <p><input type="checkbox"/> Extrinsic</p> <p><input type="checkbox"/> Intrinsic</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p>	<h3>Heel</h3> <p>Deep heel cup</p> <p style="text-align: center;"></p> <p>Heel cushion</p> <p style="text-align: center;"></p> <p>Donut cushion</p> <p style="text-align: center;"></p> <p>Medial heel skive</p> <p style="text-align: center;"></p> <p>Horseshoe cushion</p> <p style="text-align: center;"></p> <p>Medial accom.</p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p>	<h3>Metatarsal</h3> <p>Cutout</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center;">1 2 3 4 5 1 2 3 4 5</p> <p style="text-align: center; font-size: x-small;">(In cushion only)</p> <p>Met pad</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center;"><input type="checkbox"/> Soft <input type="checkbox"/> Firm</p> <p>Full met pad</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p>Met raise (in shell)*</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p>Dancer's pad</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center;">1 2 3 4 5 1 2 3 4 5</p> <p>Neuroma pad</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center;">1 2 3 4 1 2 3 4</p> <p>Shaft pad</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center;">1 2 3 4 5 1 2 3 4 5</p>	<h3>Forefoot</h3> <p>Forefoot wedge</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center;"><input type="checkbox"/> Varus <input type="checkbox"/> Valgus</p> <p>Morton's extension</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(In cushion only)</p> <p>Rev. Morton's extension</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Cutout in pad only)</p> <p>Hallux accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p>Toe crest accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p>
<p>Left</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p>	<p>Right</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p>								
<p>Left</p> <p>1/8" (3.0mm)</p> <p>2/8" (6.0mm)</p> <p>3/8" (9.0mm)</p> <p>4/8" (12.0mm)</p>	<p>Right</p> <p>1/8" (3.0mm)</p> <p>2/8" (6.0mm)</p> <p>3/8" (9.0mm)</p> <p>4/8" (12.0mm)</p>								
<p>Left</p> <p><input type="checkbox"/> Extrinsic</p> <p><input type="checkbox"/> Intrinsic</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p>	<p>Right</p> <p><input type="checkbox"/> Extrinsic</p> <p><input type="checkbox"/> Intrinsic</p> <p>Varus _____ °</p> <p>Vertical _____ °</p> <p>Valgus _____ °</p>								
<h3>Arch/midfoot</h3> <table style="width:100%; border: none;"> <tr> <td style="width:50%; border: none;"> <p>High medial overlay</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> </td> <td style="width:50%; border: none;"> <p>Base of 5th/styloid</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p> </td> </tr> <tr> <td style="border: none;"> <p>Scaphoid pad</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> </td> <td style="border: none;"> <p>Cuboid accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p> </td> </tr> <tr> <td style="border: none;"> <p>Plantar fascia accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(In cushion only)</p> </td> <td style="border: none;"> <p>Navicular accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p> </td> </tr> </table>				<p>High medial overlay</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p>	<p>Base of 5th/styloid</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p>	<p>Scaphoid pad</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p>	<p>Cuboid accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p>	<p>Plantar fascia accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(In cushion only)</p>	<p>Navicular accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p>
<p>High medial overlay</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p>	<p>Base of 5th/styloid</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p>								
<p>Scaphoid pad</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p>	<p>Cuboid accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p>								
<p>Plantar fascia accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(In cushion only)</p>	<p>Navicular accom.</p> <p style="text-align: center;"></p> <p style="text-align: center;"></p> <p style="text-align: center; font-size: x-small;">(Pad only)</p>								
 <p style="text-align: center; font-size: x-small;">Dorsal view</p>									

NOTES AND SIGNATURE

Signature: _____

Print name: _____

Date: _____